



Manuals.plus / > Hyperice / > Hyperice Hypervolt 2 Pro - Instruction Manual

Hyperice Hypervolt 2 Pro

Sign in to save this manual

Sign in to bookmark this manual and keep private notes.

Hyperice Hypervolt 2 Pro Instruction Manual

Brand: Hyperice

Advertisements

Model: Hypervolt 2 Pro

INTRODUCTION

The Hyperice Hypervolt 2 Pro is a powerful handheld percussion massage device designed to provide deep-penetrating relief to stiff muscles, aiding in faster recovery and enhanced performance. This manual provides essential information for the safe and effective use of your device.

Advertisements



WHAT'S INCLUDED

Your Hypervolt 2 Pro package includes:

- Hypervolt 2 Pro device
- 5 interchangeable head attachments: Fork, Ball, Cushion, Flat, and Bullet.
- Power supply and adapters
- Head attachment pouch

Advertisements





What's included with your Hypervolt 2 Pro

Hypervolt 2 Pro
The actual device

Power supply + adapters
To keep your device charged no matter where you are

Head attachments (5)
To tailor the device to best suit the area you're focused on

Head attachment pouch
To keep all the attachments together and organized

Figure 1: Contents of the Hypervolt 2 Pro package, including the device, various head attachments, power supply, and a storage pouch.

SETUP

Charging the Device

Before first use, fully charge the Hypervolt 2 Pro for up to four hours using only the provided Hyperice charger. Full charge is indicated when the LED light band on the handle remains green and is fully illuminated.

Advertisements



Attaching a Head

Select the desired head attachment (Fork, Ball, Cushion, Flat, or Bullet) and firmly insert it into the opening on the device. To remove, gently pull the head attachment out.



Figure 2: The Hypervolt 2 Pro device shown with its five distinct head attachments, designed for various muscle groups and intensity levels.

OPERATING INSTRUCTIONS

Powering On/Off

To turn on the device, locate the switch at the bottom of the handle and slide it to the ON position. To turn off, slide the switch to the OFF position.



Adjusting Speed Settings

The Hypervolt 2 Pro features five speed settings. Once powered on, rotate the digital speed dial located on the device to select your preferred intensity level. The LED indicators around the dial will illuminate to show the current speed setting.

Details and specifications

Weight

2.6 lbs

Size

8 x 2.5 x 10.4 in

Battery life

3 hours

Speed

5 speeds of percussion

Motor

Brushless high-torque 90W motor

Bluetooth

Connects via Bluetooth to the Hyperice App

Pressure sensor

Patented Pressure Sensor

Technology showcases how much pressure is being applied

TSA approved

Approved by TSA for carry-on



Figure 3: The digital speed dial on the Hypervolt 2 Pro, used to select from five levels of percussion intensity.

Using the Hyperice App (Bluetooth Connectivity)

For guided routines and automatic speed adjustment, connect your Hypervolt 2 Pro to the Hyperice App via Bluetooth. The app provides personalized programs to optimize your warm-up and recovery.

Advertisements



Application Tips

Apply the massager to the desired muscle group. Move the device slowly and gently, allowing the percussion to work. Avoid applying excessive pressure. The patented pressure sensor technology showcases how much pressure is being applied.



Hypervolt 2 Pro

Our most powerful device ever, now with five variable speeds and a new digital dial, the Hypervolt 2 Pro helps you recover faster so you can get back to doing more of what you love.

Figure 4: Proper application of the Hypervolt 2 Pro on the thigh for muscle recovery.

MAINTENANCE AND CARE

To ensure the longevity of your Hypervolt 2 Pro, follow these maintenance guidelines:

Advertisements



- Clean the device with a damp cloth after each use. Do not submerge the device in water.
- Store the device and its attachments in the provided pouch in a cool, dry place.
- Only use the official Hyperice charger for charging.
- Avoid dropping the device or exposing it to extreme temperatures.

TROUBLESHOOTING

Device Not Turning On

- Ensure the device is fully charged. Check the LED light band for green illumination.
- Verify the power switch at the bottom is in the ON position.

Weak Percussion/Motor Issues

- Ensure the battery is adequately charged.
- Check that the head attachment is securely fastened.
- Avoid applying excessive pressure, as the device has a built-in pressure sensor to prevent motor strain.

Advertisements



Bluetooth Connectivity Problems

- Ensure Bluetooth is enabled on your mobile device.
- Make sure the Hypervolt 2 Pro is within range of your mobile device.
- Restart both the Hypervolt 2 Pro and your mobile device.
- Ensure the Hyperice App is updated to the latest version.

SPECIFICATIONS

Feature	Detail
Weight	2.6 Pounds
Dimensions	13 x 10 x 6 inches
Battery Life	Approximately 3 hours
Speeds	5 speeds of percussion
Motor	Brushless high-torque 90W motor
Connectivity	Bluetooth (Hyperice App compatible)
Pressure Sensor	Patented Pressure Sensor Technology
TSA Approved	Yes, for carry-on
Material	Plastic

Details and specifications

Weight

2.6 lbs

Size

8 x 2.5 x 10.4 in

Battery life

3 hours

Speed

5 speeds of percussion

Motor

Brushless high-torque 90W motor

Bluetooth

Connects via Bluetooth to the Hyperice App

Pressure sensor

Patented Pressure Sensor Technology showcases how much pressure is being applied

TSA approved

Approved by TSA for carry-on



Figure 5: Visual summary of Hypervolt 2 Pro details and specifications.

SAFETY INFORMATION

Please read and understand all safety warnings and instructions before using the Hypervolt 2 Pro. Failure to do so may result in injury or device damage.

Advertisements



General Safety Precautions

- Do not operate the device in wet or damp conditions.

- Keep the device away from heat sources and direct sunlight.
- Do not insert any objects into the ventilation openings.
- This device is intended for adult use only. Keep out of reach of children.
- Consult a healthcare professional before use if you have any pre-existing medical conditions, are pregnant, or have a pacemaker.
- Do not use on injured or broken skin.
- Discontinue use immediately if you experience pain or discomfort.

Charging Safety

- Use only the charger supplied by Hyperice.
- Do not attempt to modify or repair the device or battery.
- Unplug the charger from the power outlet when not in use or before cleaning.

WARRANTY AND SUPPORT

For warranty information, product registration, and customer support, please visit the official Hyperice website or contact their customer service department. Retain your proof of purchase for warranty claims.

Advertisements



Manufacturer: Hyperice

Date First Available: October 12, 2021

Advertisements



Related Documents - Hypervolt 2 Pro

	Hypervolt Go 2 Percussion Massage Device: Operating Instructions and Safety Guide Comprehensive guide to the Hypervolt Go 2 handheld percussion massage device by Hyperice, covering operating instructions, safety warnings, charging, cleaning, specifications, and warranty information.
	Hypervolt Operating Instructions and Safety Guide Official operating instructions and safety guidelines for the Hypervolt percussion massager by Hyperice. Learn how to use, clean, and maintain your device safely.

	Hypervolt BT Operating Instructions and Safety Information Comprehensive operating instructions, safety guidelines, and specifications for the Hypervolt BT percussion massage device. Learn how to use, clean, and store your Hypervolt BT for optimal muscle recovery and circulation.
	Hypervolt GO Operating Instructions and Safety Guide Comprehensive operating instructions, safety warnings, specifications, and usage guidelines for the Hypervolt GO handheld percussion massage device by Hyperice.
	Hypervolt GO User Manual and Instructions Discover the Hypervolt GO, a handheld percussion massage device designed to relieve muscle soreness and stiffness, accelerating warm-up and recovery post-activity. This manual provides operating instructions, safety guidelines, specifications, and care information.
	Hypervolt 2 Percussion Massage Device: Operating Instructions and Safety Manual Get detailed operating instructions and safety guidelines for the Hypervolt 2 percussion massage device by Hyperice. Learn how to use, charge, and maintain your device safely.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

HYPERICE HYPERVOLT 2 PRO

KEEP EXPLORING MANUALS, GUIDES, AND SUPPORT DOCS.

Search model numbers, browse help resources, and jump straight into the next PDF manual or troubleshooting page without losing your place.

Hyperice Hypervolt 2 Pro

Search Manuals+

EXPLORE

Deep Search

Upload

Login

MANUALS+

Home

About Us

Privacy Policy

Login

GUIDES

Amazon Manuals

AliExpress Manuals

Video Guides

Login

FOLLOW

@manualsplus

YouTube

LinkedIn

Facebook

Login

Manuals+

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners.

SAVE MANUALS, NOTES, UPLOADS, AND QUESTIONS

Sign in with Google, X, or a passkey to keep private notes on saved manuals, track public uploads, and revisit questions you posted.

Open login options

[Privacy Policy](#) · [DMCA](#)